

<b>IHE Integration Statement</b>		<b>Date: 2005-05-30</b>
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>
CHILI GmbH	CHILI/PACS	3.5
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Consistent Time	Time Client	None
Consistent Time	Time Server	None
Patient Information Reconciliation	Image Manager/Archive	None
Patient Information Reconciliation	Performed Procedure Step Manager	None
Basic Security	Image Manager/Archive	None
Scheduled Workflow	Image Manager/Archive	None
Scheduled Workflow	Performed Procedure Step Manager	None
<b>Internet address for vendor's IHE information :</b> <a href="http://www.chili-radiology.com/ihe">http://www.chili-radiology.com/ihe</a>		
<b>Links to Standards Conformance Statements for the implementation</b>		
<b>Links to general information on IHE</b>		
<b>In North America :</b> <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>	<b>In Europe :</b> <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	<b>In Japan :</b> <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>