

IHE Integration Statement		Date: 2006
Vendor	Product Name	Version
CHILI GmbH	CHILI PACS	3.6
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Report Repository	None
Access to Radiology Information	External Report Repository Access	None
Basic Security	Image Manager/Archive	None
Scheduled Workflow	Department System Scheduler/Order Filler	None
Scheduled Workflow	Image Display	None
Scheduled Workflow	Image Manager/Archive	None
Patient Information Reconciliation	Department System Scheduler/Order Filler	None
Patient Information Reconciliation	Image Manager/Archive	None
Consistent Presentation of Images	Print Composer	None
Portable Data for Imaging	Image Display	None
Portable Data for Imaging	Print Composer	None
Portable Data for Imaging	Portable Media Creator	None
Portable Data for Imaging	Portable Media Importer	None
Consistent Time	Time Client	None
Consistent Time	Time Server	None
Internet address for vendor's IHE information : http://www.chili-radiology.com/ihe		
Links to Standards Conformance Statements for the implementation		
Links to general information on IHE		
In North America : http://www.ihe.net	In Europe : http://www.ihe-europe.org	In Japan : http://www.jira-net.or.jp/ihe-j