

IHE Integration Statement		Date: 2005-05-30
Vendor	Product Name	Version
CHILI GmbH	CHILI/Report	3.5
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Consistent Presentation of Images	Print Composer	None
Consistent Time	Time Client	None
Scheduled Workflow	Image Display	None
Scheduled Workflow	Image Manager/Archive	None
Scheduled Workflow	Performed Procedure Step Manager	None
Internet address for vendor's IHE information : http://		
Links to Standards Conformance Statements for the implementation		
Links to general information on IHE		
In North America : http://www.rsna.org/IHE	In Europe : http://www.ihe-europe.org	In Japan : http://www.jira-net.or.jp/ihe-j