
CHILI PACS

IHE Integration Statement Version 2019



Manufacturer CHILI GmbH

Friedrich-Ebert-Str. 2
D-69221 Dossenheim/Heidelberg
Tel. (+49) 6221 - 1 80 79 -10
Fax (+49) 6221 - 1 80 79 -11

Support: (+49) 6221 - 1 80 79 -20
E-mail: support@chili-radiology.com
Internet: www.chili-radiology.com

Revised: April 2019
CHILI PACS Version 4.64
IHE Integration Statement Version 2019



© 2003–2019 CHILI GmbH

The authors reserve the right to make changes at any time without prior notice.
All rights reserved, in particular the right to reproduction, distribution, and translation. No part of the software or the user manual may be reproduced or processed, copied, or distributed by electronic means without the prior permission in writing from the authors.

www.chili-radiology.com



CHILI®

CHILI AND THE CHILI-LOGO ARE REGISTERED
TRADEMARKS OF THE CHILI GMBH.

| IHE Integration Statement | | Date: 2019-04-17 |
|--|-----------------------------------|---|
| Vendor | Product Name | Version |
| CHILI GmbH | CHILI PACS | 4.64 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integrations Profiles, Actors, and Options listed below: | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented |
| Access to Radiology Information (ARI) | Report Repository | None |
| | External Report Repository Access | None |
| Audit Trail and Node Authentication (ATNA) | Secure Node | None |
| | Secure Node | Radiology Audit Trail |
| | Secure Node | FQDN Validation of Server Certificate |
| | Secure Node | BCP195 TLS Secure Transport Connection - All TLS Versions |
| | Secure Node | BCP195 TLS Secure Transport Connection - TLS 1.2 Floor |
| Basic Security (SEC) | Image Manager/Archive | None |
| Consistent Presentation of Images (CPI) | Print Composer | None |
| | Image Manager | None |
| | Image Display | None |
| Consistent Time (CT) | Time Client | None |
| | Time Server | None |
| Cross-Enterprise Document Media Interchange (XDM) | Portable Media Creator | CD-R, USB |
| | Portable Media Importer | CD-R, USB |
| Cross-Enterprise Document Sharing for Imaging (XDS-I.b) | Imaging Document Consumer | None |
| | Imaging Document Source | CDA Wrapped Text Report |
| | Imaging Document Source | PDF Report |
| | Imaging Document Source | Set of DICOM Instances |
| Evidence Documents (ED) | Image Display | None |
| | Image Manager/Archive | None |
| Healthcare Provider Directory (HPD) | Provider Information Directory | Provider Information Feed |
| Imaging Object Change Management (IOCM) | Image Display | None |
| | Image Manager/Archive | None |
| | Change Requester | None |
| Import Reconciliation Workflow (IRWF) | Importer | Scheduled Import - uses MWL, Unscheduled Import - no MWL |
| | Image Manager | None |
| Import Reconciliation Workflow (IRWF.b) | Importer | None |
| | Image Manager/Archive | None |
| Invoke Image Display (IID) | Image Display | None |
| Key Image Note (KIN) | Image Manager/Archive | None |
| | Evidence Creator | None |
| | Image Display | None |
| Mammography Image (MAMMO) | Image Manager | None |

| Integration Profiles Implemented | Actors Implemented | Options Implemented |
|--|---|---------------------|
| Patient Administration Management (PAM) | Patient Demographics Consumer | None |
| | Patient Demographics Consumer | Merge |
| Patient Demographic Query (PDQ) | Patient Demographics Consumer | None |
| Patient Identifier Cross-referencing for MPI (PIX) | Patient Identifier Cross-reference Consumer | None |
| Patient Information Reconciliation (PIR) | Department System Scheduler/Order Filler | None |
| | Image Manager/Archive | None |
| Portable Data for Imaging (PDI) | Image Display | None |
| | Print Composer | None |
| | Portable Media Creator | None |
| | Portable Media Importer | None |
| Radiation Exposure Monitoring (REM) | Image Display | None |
| | Image Manager | None |
| Scheduled Workflow (SWF) | Department System Scheduler/Order Filler | None |
| | Image Display | None |
| | Image Manager/Archive | None |
| | Performed Procedure Step Manager | None |
| Scheduled Workflow (SWF.b) | Department System Scheduler/Order Filler | None |
| | Image Display | None |
| | Image Manager/Archive | None |
| | Performed Procedure Step Manager | None |
| Simple Image and Numeric Report (SINR) | Report Manager | None |
| Teaching File and Clinical Trial Export (TCE) | Export Manager | None |
| | Export Selector | None |
| Web-based Image Capture (WIC) | Image Capturer | JPEG Storage |
| | Image Manager | None |
| Web-based Image Access (WIA) | Imaging Document Responder | None |
| | Imaging Document Source | None |

| | | |
|---|--|---|
| Internet address for vendor's IHE information: http://www.chili-radiology.com/ihe | | |
| Links to Standards Conformance Statements for the Implementation: http://www.chili-radiology.com/ihe/ | | |
| Links to general information on IHE | | |
| In North America: http://www.ihe.net | In Europe: http://www.ihe-europe.net | In Japan: http://www.ihe-j.org |